
DOC022.97.80411

# QuikChem<sup>®</sup> 8500 Series 2

01/2016, Edition 2

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Troubleshooting on page 14

## General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

## Use of hazard information

### ▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### ▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### ▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the need for protective eye wear.
	This symbol indicates that the marked item can be hot and should not be touched without care.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.
	This symbol indicates that the object is heavy.

## Omnion software installed

This instrument is used with the Omnion software installed on the computer. For specific information, refer to the Omnion documentation.

## Operation

### Operation overview

<b>▲ WARNING</b>	
 	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

## ▲ WARNING

Fire hazard. The user is responsible to make sure that sufficient precautions are taken when the equipment is used with methods that use flammable liquids. Make sure to obey correct user precautions and safety protocols. This includes, but is not limited to, spill and leak controls, proper ventilation, no unattended use, and that the instrument is never left unattended while power is applied.


The QuikChem® System automatically does wet chemical determinations with the principle of Flow Injection Analysis. The pump pulls the sample from the sampler and to the injection valve. At the same time, reagents are continuously pumped through the system. The sample is loaded in the sample loop of one or more injection valves. The injection valve then connects the sample loop in line with the carrier stream. This moves the sample out of the sample loop and on the manifold.

The sample and reagents then come together in the manifold (reaction module) where the sample can be diluted, concentrated, dialyzed, extracted, incubated and derivatized. Mixing occurs in the narrow bore tubing under laminar flow conditions. For each method, the operating parameters are optimized to give high sample throughput, high precision and high accuracy. The FIA peaks are transient and show the non-steady state conditions normally used. In FIA, there is usually no equilibrium.

### Manifold diagram and symbols

Figure 1 shows an example of a manifold diagram for the analyzer. The manifold diagram can show a back pressure loop. Refer to the maintenance and troubleshooting manual for more information.

Figure 1 Manifold diagram


1 Sample (green)	7 Reagent connection to manifold
2 Carrier solution	8 Probe rinse
3 Pump flow	9 Flow cell
4 Probe wash	10 To waste
5 Reagent	11 To port 6 of the next valve or waste
6 Reagent	12 Pump

### Mixing coil symbol

The number next to the mixing coil symbol shows the length of the coil support. The length of the coil support is measured from hole to hole. Mixing coils typically have a 0.8 mm ID Teflon tubing that winds around a coil support. The length of Teflon tubing to wind on the coil support is given in each QuikChem method. Refer to Figure 2.


**Figure 2 Mixing coil symbol**

### Injection valve symbols

Refer to [Figure 3](#). The injection valve is on the left side of every channel of the system unit. The numbers around the hexagon show the six ports on the valve. The zig-zag or the coil symbols show the sample loops. The sample loop tubing is connected between ports 1 and 4. The length of this tubing is shown on the manifold diagram. Refer to the QuikChem method. The lines that come out of the hexagon are the Teflon tubing connected to the valve ports.

**Figure 3 Injection valve symbols**


1 Primary

2 Secondary

### Pump and pump flow symbols

[Figure 4](#) shows the reagent pump. The top arrow shows the pump flow direction. The line below this refers to the rinse water line attached to the wash reservoir in the sampler. The wash line is connected to the bottom nipple of the wash reservoir.

**Figure 4 Pump and pump flow symbols**


1 Sample (green or purple)	4 To wash bath fill
2 Carrier	5 Color reagent
3 From water pump tube	6 Pump tube (color)

### Reagent line representation

The reagents and sample lines are shown by a straight line that goes across the rectangle. Refer to [Figure 5](#). The specified color is the color of the pump tube tabs to use for that line. All pump tubes have different color tabs on each side to show the inner diameter of that tubing. The reagent tubing has a label that shows the solution that goes through them. The carrier tubing is always connected to port 2 of the valve. The sample tubing is always connected to port 6 of the valve. The sample tubing is always connected to port 2 of the valve.

**Figure 5 Reagent line representation**


1 Sample (green)	3 Pump tube (color)
2 Carrier	4 Color reagent

### Tee connector symbol

Every intersection of two reagent tubes shows a tee connector. Refer to [Figure 6](#). Tee connectors mix two reagents. The tee has three ports. Two side ports (inlets) are connected to two reagent tubes. The top port (outlet) can be connected to a mixing coil or an item shown in the manifold diagram.

**Figure 6 Tee connectors symbols**


1 Tee connector fitting	3 Alternate symbol
2 Tee connector symbol	

**Heater block symbol**

Figure 7 shows a heater block installed on the system unit. Set the heater at the temperature specified in the diagram (60 °C (140 °F) in the example). The QuikChem method shows the length of tubing to wind on the heater block.


**Figure 7 Heater block symbols**

**Injection valve operation**

The valve has a stator and a rotor. Refer to Figure 8. The rotor has three pathways. The liquid moves through the pathways. The alignment of the pathways with the ports on the stator control the direction in which the liquid flows through the valve. Refer to Table 1 for the ports on the six-port injection valve. The injection valve has two positions: the inject position and the load position. Refer to Inject position on page 9 and Load position on page 9. The valve moves bi-directionally between those two positions. When the system is not in operation, the valves stay in the inject position.

**Figure 8 Stator and rotor**

1 Stator	2 Rotor
----------	---------

**Table 1 Port connections**


Port connection	Description
Sample tubing	The sample moves from the sampler or the sample probe to port 6 on the injection valve.
Sample loop	There is tubing between ports 1 and 4. The sample loop finds the volume of sample to be put on the manifold for analysis.
Carrier tubing	The carrier is an artificial blank solution that will move the sample on the manifold. The carrier tubing input is always connected to port 2. The carrier solution and the carrier tubing to be used for an analysis are specified in the QuikChem method.
Manifold	The carrier stream exits the valve at port 3 and flows to the manifold.
Waste tubing	When only one channel operates, the waste tubing is connected to port 5 on the injection valve. When two or more channels operate simultaneously, port 5 of the first injection valve is connected to port 6 of the second injection valve. Port 5 of the second injection valve is connected to waste or to the next valve in the series. The connection between ports 5 and 6 have 15 cm Teflon® tubing 0.8 mm ID.

### Inject position

At the inject position, the carrier flows through the sample loop to the manifold. When the valve moves to the load position, some of the carrier is caught in the sample loop. The carrier that is caught in the sample loop at the inject position is moved out to waste by the sample stream. The sample stream moves through the sample loop and goes to waste.

When the valve is back in the inject position, the carrier stream flows to port 2, then to port 1 and through the sample loop. The carrier moves the sample out through port 3 and to the manifold. The sample stream goes in port 6 and out of port 5 to waste. Refer to [Figure 9](#).

**Figure 9 Inject position**


1 Waste tubing/next valve	4 Carrier line in
2 Sample in from sampler or adjacent valve	5 Flow to manifold
3 Sample loop	

### Load position

At the load position, the sample loop is filled with sample. When the valve moves back to the inject position, the sample that flows through the sample loop is caught momentarily.

The sample volume can be calculated from the length and the inner diameter of the tubing used. The sample is moved on the manifold by carrier for reaction and analysis. Refer to [Figure 10](#).

**Figure 10 Load position**

1 Waste tubing/next valve	4 Carrier tubing
2 Sample in from sampler or adjacent valve	5 Flow to manifold
3 Sample loop	

## Start up the system unit

Connect all the modules of the system to a power strip. Do not set the pump to on until the instrument is ready to operate.

1. Set the power strip to on. Then set each individual module switch to on.  
**Note:** A 6-second delay occurs from the power circuits for the instrument to start up. At this time, a beeping sound is audible.  
The sampler automatically completes an operation check. The probe moves over the wash bath or directly above the wash reservoir. When the sampler is initialized through the sampler configuration screen or when a run is started, the probe goes in the wash reservoir.
2. Look through the side panel on the system unit to see if the lamp is on. The lamp comes on in about one minute. The delay lengthens for the life span of the lamp. The heater controllers shows the current temperature of the blocks.
3. Start the software to make sure that the software is connected to the system unit. A dialog box shows 'Connecting to System unit(s)'. If the setup was done, the dialog box goes away.  
**Note:** If the connection is not made correctly, a dialog box shows a message that system unit A is not found. Make the connections and click **OK**.

## Prepare the instrument for analysis

When the system start up is done, prepare the instrument for analysis. Refer to [Start up the system unit](#) on page 10.

1. Make sure to use a special detector for the analysis. A photometric detector module is necessary for most QuikChem methods.  
**Note:** The Ion Chromatography (IC) channel (optional) typically uses a conductivity module instead of a photometric detector module.
2. Install the manifold on the channel that will do the analysis. For the manifold installation, refer to the installation manual.
3. Make all the injection valve fluidic connections. Refer to [Injection valve operation](#) on page 8.

4. Make all the flow cell fluidic connections to the detector.
5. Set all the pump tubes on the pump. The number of pump tubes varies from method to method, but the sample tube and wash tube (for the sampler) must be used for analysis.
6. Flush deionized water through all the tubes to make sure that there are no leaks. When a column is necessary (e.g., for a nitrate analysis or other method), make sure that the column is offline. Remove the flow cell from the detector and look for leaks. If no leaks are found, put the reagents online.  
*Note: Make sure to read the QuikChem method notes for the startup of specific methods.*
7. Set the heater controller temperature in the Omnion software under the Analyte tab. Refer to the QuikChem method for temperature settings. Refer to the online help for heater configuration information.
8. If the system has a dilutor, install the diluent line E in the diluent (i.e., deionized water). Put the empty test tubes in the empty tubes rack.
9. Pour the calibration standards in the standard vials. Pour some samples in the test tubes.
10. Open a 'run worksheet' or open a new worksheet. The information in the worksheet (standards and samples) tells the instrument what analysis to do and starts the analysis. Refer to [Calibrate](#) on page 13.

## Operate the pump

To use new pump tubing for the first time, operate the pump 10 minutes to get constant flow rates. [Table 2](#) shows the functions of each button on the pump.

1. Put all the reagent tubing in the corresponding containers or deionized water.
2. Before a reagent is used, flush deionized water through all of the lines to find leaks.
3. If the manifold has a cadmium column, make sure that it is bypassed. Water and air will cause damage to the cadmium column.  
*Note: For some methods, it is necessary to move some reagents first. For example, the TKN method has a buffer solution that must flow first to prevent a white precipitate in the manifold tubing. Refer to the QuikChem method for instructions.*

**Table 2 RP-150 pump button functions**

Button	Function
MIN	Sets the pump speed to 4.
MAX	Sets the speed to a maximum of 999. The button must be kept pushed to keep the maximum speed. <i>Note: If this button is pushed too long, the pump will stop to prevent motor burn out.</i>
NORMAL RUN	Normal operation starts the pump and lets the remote control of the pump standby speed as given by the Omnion software at the end of an operation or batch. The pump operates at the shown speed. Normal speed is 35. Use the <b>UP</b> and <b>DOWN</b> arrows to change this setting. If the Omnion software is set to Standby speed, the pump operates very slowly and the display shows 'rc 3' (remote control speed 3).
MANUAL RUN/STOP	Manual operation lets the pump override the remote control of the standby speed. This button toggles between manual speed control and manual stop. If the pump operates, this button always stops the pump. If the pump does not operate, this button still overrides the remote control from the computer and runs the pump at the set speed. The MANUAL RUN button is commonly used to make the pump go to the normal speed to fill the tubing with reagents.
UP arrow	Increases the pump speed setting. This button is not on in the MIN or MAX modes.
DOWN arrow	Decreases the pump speed setting. This button is not on in the MIN or MAX modes.

## Shut down the pump

Do the steps that follow to shut down the system after an operation and/or end of the day.

1. Rinse the reagent lines with deionized water.
2. Push the **MANUAL RUN/STOP** button to stop the pump.
3. Set the green power switch or the power strip to off.
4. Release all the cartridges. If the pump tubes pushed down for more than a few minutes, it causes the tubing to become flat and must be replaced.
5. Push on the cartridge holder until the pump tube cartridge is released from the pump. The pump tube can stay on the cartridges. The tension lever has no effect on the pump tube when the cartridge is not pushed down.

## Prepare a universal dye

To examine the operation of the system, operate the tests with the universal dye instead of reagents. Make sure that the universal dye that comes with the system has been diluted with the directions on the label.

**Note:** *The dye will not be found with the 880 nm filter.*

## Temperature controller

Refer to the Omnion software to set the heater properties with the operating temperature (set-point), when to send the set-point to the heater element and post-run set-point. The temperature can be seen real-time in the software and from the controller on the heater assembly. Refer to the Omnion help file for more information.

**Note:** *Always set the run and the post-run temperature in the Omnion software.*

## Prepare the dilutor

For systems that use the PDS200 dilutor, make sure to put the diluent line in the diluent (usually deionized water). All triggered dilutions will be completed at the end of the tray run. All necessary dilutions will be completed at the time the sample is in the tray run sequence.

**Triggered dilutions:** The instrument recognizes a sample that is over the concentration range of the method. The default limit is 10% over the range of the highest concentration standard. Above this limit, the dilutor automatically dilutes all over-range samples. This feature must be selected in the analyte table in the Omnion software.

**Requested dilutions:** Set the dilution factor for a specific sample. This factor is specified by the Run Worksheet. Refer to the Omnion help guide for more information.

**Note:** *Requested dilutions can be used to make calibration standards from a stock standard. Make sure that the tray includes all the calibration standards. Use the same cup number for all of them (e.g., cup 15). Put the stock standard in that cup (e.g., 15). Show the dilution factor for each calibration standard accordingly.*

The dilutor system uses an empty tubes rack to complete dilutions. The ASX-520 series sampler has the fourth rack that is reserved for empty tubes.

## Set the valve timing

Make sure that the valve timing is set correctly. Valve timing that is not correct can result in the absence of peaks, presence of air spikes or very small inaccurate peaks. Set the cycle period parameter as specified in the QuikChem method. When multiple methods are operated simultaneously, use the longest cycle period in the QuikChem method. Do the steps that follow to make sure that the valve timing is correct.

**Note:** *The "Time to Valve" value of channels 2 and higher is measured from when the probe that reaches the dye, through all previous sample loops until it arrives at port 6 of the channel. A "Time to Valve" value must be entered for all channels that are used.*

1. Put in the universal dye as the sample.
2. Look at the sample slug and record the amount of time for it to get through the test tube and to port 6. Record the time when the sample probe goes in the test tube.
3. Record the time when the dye reaches port 6 at the valve that is measured. This time is called "Time to Valve".

## Set the sampler timing

When possible, set the sampler timing at the values specified in the Lachat method. Do the steps that follow when less sample is necessary for the analysis.

1. Decrease the sample period.
2. Examine the new Run Timing setting with dye. Refer to [Prepare a universal dye](#) on page 12. If either of these conditions occur, the computer will show an error message:
  - The Sample Period plus Minimum Probe in Wash Period is longer than the Method Cycle Period
  - The Load plus Injection Period is longer than the Method Cycle Period

## Calibrate

After the manifold installation, refer to the QuikChem method for the calibration details.

1. Set the power to on. Make sure that all system components are set to on.
2. Open the Omnion software.
3. Open a run file.
4. Prepare the reagents.
5. Install the manifold.
6. Flush deionized water through the manifold.
7. Look for leaks on the system. If no leaks are found, operate the system with reagents.
8. Prepare the calibration standards.
9. Install the calibration standards and the samples on the sampler.
10. Prepare the concentration in Run properties>Samples.
11. Start the run with the calibration standards and samples.
12. After the analysis is complete, print/export the report. For print options, refer to the online help.
13. Start the shutdown procedure to rinse the system. Refer to [Shut down the system](#) on page 13.

## Shut down the system

1. To use a manifold with a column, such as the cadmium column for nitrate, make the column offline. Set the switching valve to the bypass position. If the manifold has a column, but does not have a switching valve, set the pump power to off and replace the column with a piece of Teflon tubing. Start the pump.  
*Note: Some methods make it necessary for some reagents (typically, the buffer) to be removed last. Look at the QuikChem method for instructions.*
2. Remove the reagent tubing from each reagent and rinse the tubing and glass weights. Keep the lines and weights in rinse solutions.
3. If the QuikChem method recommends a rinse solution, install all reagent transmission tubing in the solution and operate the pump for five minutes at standard speed.  
*Note: This is a critical step in the preventive maintenance of the manifolds.*
4. Install the transmission tubing in deionized water and let the system rinse for 10 minutes at standard speed.
5. If the manifold will not be used for more than three days or it is removed and kept in storage:
  - a. Remove the transmission tubing from the deionized water.
  - b. Use the pump to remove all liquid from the manifold.
6. Set the power to the pump to off. Release the pump tube cartridges.
7. Close all files on the computer.

- Set the power to the power strip to off.

**Note:** If a lower set-point temperature is not given in the software, the heater controller stays at the original set temperature. If the method makes it necessary for a temperature higher than 60 °C (140 °F), lower the temperature setting to 60 °C (140 °F) or less before power is removed from the system.


## Troubleshooting

### Do a pump test

Do the steps that follow to make sure that the pump speed is correct.

- Use an accurate watch or chronometer to measure the time for the pump to make 10 revolutions.
- Remove all pump cartridges from the pump. Look for the index marks as shown in [Figure 11](#). Index marks let the pump revolutions be counted. Manually move the rollers to find the index mark on the roller mounting. If either mark is hard to read, make new index marks with a marker.
- Attach a new green-green pump tubing. Make sure that the tension lever is set correctly at the 12 o'clock (vertical) position.
- Set the pump power to on. Push **NORMAL RUN**. The display must show a speed of 35 to continue this test. If the display shows anything else, refer to the installation manual.
- Use an accurate watch or chronometer to measure the time for the pump to make 10 revolutions. Make sure the time is 50 seconds  $\pm$  1 second.
- If the recorded time is less than 49 seconds or more than 51 seconds, contact technical support or a local sales representative for assistance.

**Figure 11** Index marks on the pump

1 Index mark that turns with the rollers

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)