



ISO 9001 Certified

TEST CERTIFICATE

Model *Pocket Colorimeter II*

Part Number 58700-00

Serial Number 14010E239165

Inspected By Pocket Colorimeter Team

Date 01/14/2014

CHARACTERISTIC	STANDARD	ACTUAL
Visual Inspection	PASS/FAIL	PASS
Programming	PASS/FAIL	PASS
Display Test	PASS/FAIL	PASS
Photometric Accuracy	1.050 ±0.100 Abs	1.072

*The photometric accuracy specification includes uncertainties due to filter values, reproducibility of filter insertion and instrument readings, and photometric linearity error.

Test equipment used for the verification of Hach manufactured instruments is calibrated using standards traceable to National Institute of Standards and Technology (NIST) standards. Where such standards do not exist, the basis for calibration is documented.

Production Engineering Approved By Steve Kemper

Date 03/11/2014

Customer Name: _____

Customer P.O.# _____

Location: _____

Hach Ref. # _____

THIS IS AN EXAMPLE CERTIFICATE